

**Product range:**

- MiniCollect® Tube FX Sodium Fluoride / Potassium Oxalate
- MiniCollect® Complete FX Sodium Fluoride / Potassium Oxalate

**General description of the product:**

The tubes contain potassium oxalate and sodium fluoride.

The cap colour is grey and the tubes are used for capillary blood.

**Field of application:**

MiniCollect® FX Sodium Fluoride/ Potassium Oxalate Tubes are used for the determination of glucose and lactate in capillary blood.

**Intended purpose (Use)**

MiniCollect® FX Sodium Fluoride/ Potassium Oxalate Tubes are used for the determination of glucose and lactate in capillary blood.

**Intended user:**

MiniCollect® Tubes have to be used by qualified personnel only.

**Groups of patients:**

No restrictions.

**Device configurations and variations**

In Table the variations of MiniCollect® FX Sodium Fluoride / Potassium Oxalate Tubes and MiniCollect® Complete FX Sodium Fluoride / Potassium Oxalate Tubes are listed.

Item Number	Product name / description
<b>450540</b>	MiniCollect® TUBE 0.25 ml FX Sodium Fluoride / Potassium Oxalate grey cap
<b>450541</b>	MiniCollect® TUBE 0.5 ml FX Sodium Fluoride / Potassium Oxalate grey cap
<b>450552</b>	MiniCollect® Complete 0.25 ml FX Sodium Fluoride / Potassium Oxalate grey cap, pre-assembled with Carrier Tube 13x75

**Shelf life**

The following table lists the shelf life of the products at a storage temperature between 4°C and 25°C.

Item number	Product	Shelf life
<b>450540</b>	MiniCollect® FX Sodium Fluoride / Potassium Oxalate – 0,25ml	534 days
<b>450541</b>	MiniCollect® FX Sodium Fluoride / Potassium Oxalate – 0,5ml	534 days
<b>450552</b>	MiniCollect® Complete FX Sodium Fluoride / Potassium Oxalate – 0,25ml	534 days

**Centrifugation**

--	--	--

<b>Tube Type</b>	<b>Inversions (mixing)</b>	<b>Recommended g-force relative centrifugal force (rcf)</b>	<b>Time (min)</b>
MiniCollect® FX Sodium Fluoride / Potassium Oxalate	8-10x	3000 g	10
MiniCollect® Complete FX Sodium Fluoride / Potassium Oxalate	8-10x	3000 g	10

Centrifugation should be done at a temperature of 15°C – 24°C.