

General Information

MiniCollect® Tubes are plastic, non-evacuated, non-sterile low sample volume tubes with integrated collection devices. The closure is color coded to identify the additives which are present in varying concentrations depending on the tube type and stated volumes.

Components

The tube consists of the following components:

- MiniCollect® Tube: With an integrated scoop for easier blood collection.
- Cap with pierceable membrane
- Additive
- MiniCollect® Carrier Tube
- Label with filling mark(s).

| Component | Material |
|---|-----------|
| MiniCollect® Cap | PE TPE |
| MiniCollect® Tube | PP |
| MiniCollect® Carrier Tube 13/75mm (only for MiniCollect® Complete) | PET |
| Tube label | PP |

Cautions / Precautions

- Do not use tubes if foreign matter is present.
- Handle all biological samples and blood collection devices according to the policies and procedures of your facility.
- Obtain appropriate medical attention in the case of any exposure to biological samples.

- Check all tubes to verify appropriate product and shelf-life before use. Do not use tubes after the expiration date.
- It is the laboratory's ultimate responsibility to verify that a change from one tube to another does not significantly affect analytical results obtained from patient samples.
- Single use only.

Sterility

Not applicable, since the product is not sterile.

Labelling

| Symbol | Title of symbol | Norms / Directives | Tube label | Rack label | Inner carton label | Outer carton label |
|-------------------------------------|-----------------------|--------------------------|------------|------------|--------------------|--------------------|
| <input checked="" type="checkbox"/> | Catalogue number | EN 980 ISO 15223-1 | X | X | X | X |
| <input checked="" type="checkbox"/> | Batch code | EN 980 ISO 15223-1 | X | X | X | X |
| <input checked="" type="checkbox"/> | Use-by date | EN 980 ISO 15223-1 | X | X | X | X |
| <input checked="" type="checkbox"/> | Communauté Européenne | Directive 98/79/EC | X | X | X | X |
| | Do not re-use | EN 980 ISO 15223-1 | X | X | X | X |

| | | | | | | |
|-------------------------------------|------------------------------------|-----------------------------|---|---|---|---|
| <input checked="" type="checkbox"/> | In vitro diagnostic medical device | EN 980 ISO 15223-1 | X | X | X | X |
| | Temperature limitation | EN 980 EN ISO 15223-1 | - | X | X | X |
| <input checked="" type="checkbox"/> | Consult instructions for use | EN 980 EN ISO 15223-1 | - | X | - | - |
| <input checked="" type="checkbox"/> | Manufacturer | EN 980 EN ISO 15223-1 | - | X | X | X |

UDI attributes

For the UDI-marking the GS1 UDI barcode is used. A GS1-128 barcode is printed on the carton label and a GS1 DataMatrix is printed on the rack label.

Both codes consist of:

(01) GTIN of the packaging unit, (17) expiration date, (10) batch number

1. GTIN (Global Trade Item Number) of the packaging unit – 14 digits (01)
2. Expiration date – 6 digits (YYMMDD) (17)
3. Batch number – 8 characters (10)

Labeling the product with a UDI Barcode is only required by the US-FDA so far.

Packaging

| | | | |
|--------------|-----------------|--------------------------|-----------------|
| MiniCollect® | | MiniCollect® Complete | |
| Packaging | Amount of tubes | Packaging | Amount of tubes |
| | | | |

| | | | |
|------------------|------------------------------|------------------|------------------------------|
| Rack | 50 Tubes | Rack | 50 Tubes |
| Inner box | 500 Tubes (= 10 racks) | Inner box | 500 Tubes (= 10 racks) |
| Outer box | 2000 Tubes (= 4 inner boxes) | Outer box | 2000 Tubes (= 4 inner boxes) |

Shipping stability - DROP TEST

The function of the package was verified in a drop test under transport conditions. Results: The outer boxes and the inner boxes did not tear or break and the product was not damaged. The packaging has passed all requirements of the drop test.

Recommended Order of Draw (based on CLSI GP42-A6)

If multiple specimens are to be collected, including EDTA specimens, the EDTA specimen is drawn first to ensure adequate volume and accurate haematology test results. Specimens with other additives are collected next; serum specimens are collected last.

1. EDTA
2. Heparin / Heparin Sep
3. Glycolytic inhibitor tubes
4. Serum / Serum Sep

Literature

GP42-A6 Procedures and Devices for the Collection of Diagnostic Capillary Blood Specimens; Approved Standard-Sixth Edition.

