

Kremsmünster, 05.02.2020

## VACUETTE® Urine Collection System – Product Data Specifications

<b>Item N°<sup>1</sup></b>	454141
<b>Characterisation<sup>1,2</sup></b>	VACUETTE® TUBE 3 ml Z Urine No Additive
<b>Description<sup>1</sup></b>	<i>Dimension</i> 13x75mm <i>Capacity</i> 3ml
<b>Materials<sup>1,3</sup></b>	<i>Tube</i> PET
<b>Additive<sup>1,2</sup></b>	none
<b>Manufacturer<sup>2</sup></b>	Greiner Bio-One GmbH Bad Haller Str. 32 A-4550 Kremsmünster
<b>Sterilisation<sup>1</sup></b>	Sterile SAL 10 <sup>-6</sup>
<b>Directive<sup>2</sup></b>	IVDD 98/79/EC “Directive 98/79/EC of the European Parliament and council of 27th of October 1998 on in vitro diagnostic medical devices“. Classification: Other device (all devices except Annex II and self-testing devices).
<b>Storage<sup>1</sup></b>	4-25°C
<b>Shelf Life<sup>2</sup></b>	none
<b>Packaging<sup>1</sup></b>	24 racks of 50 pcs., 1200 pcs. in total

Sources:

1 AT\_VAC\_LMDS\_454141\_rev00.pdf  
2 SAP



Spec 454141\_R00 page 1/2

Printed Copies will not systematically be updated. Valid at the date of issue. Information is subject to change.

Kremsmünster, 05.02.2020

**Label information<sup>1</sup>**

**Rack – Carton - Tube – Examples**



Sources:  
 1 AT\_VAC\_LMDS\_454141\_rev00.pdf  
 2 SAP

Spec 454141\_R00 page 2/2

Printed Copies will not systematically be updated. Valid at the date of issue. Information is subject to change.