

Kremsmünster, 05.02.2020

VACUETTE® Urine Collection System – Product Data Specifications

Item N°¹	456065
Characterisation^{1,2}	VACUETTE® TUBE 6 ml Z Urine No Additive
Description¹	<i>Dimension</i> 13x100mm <i>Capacity</i> 6ml
Materials^{1,3}	<i>Tube</i> PET
Additive^{1,2}	none
Manufacturer²	Greiner Bio-One GmbH Bad Haller Str. 32 A-4550 Kremsmünster
Sterilisation¹	Sterile SAL 10 ⁻⁶
Directive²	IVDD 98/79/EC “Directive 98/79/EC of the European Parliament and council of 27th of October 1998 on in vitro diagnostic medical devices“. Classification: Other device (all devices except Annex II and self-testing devices).
Storage¹	4-25°C
Shelf Life²	none
Packaging¹	24 racks of 50 pcs., 1200 pcs. in total

Sources:

1 AT_VAC_LMDS_456065_rev10.pdf
2 SAP




Spec 456065_R00 page 1/2

Printed Copies will not systematically be updated. Valid at the date of issue. Information is subject to change.

Kremsmünster, 05.02.2020

Label information¹

Rack – Carton - Tube – Examples



Sources:
 1 AT_VAC_LMDS_456065_rev10.pdf
 2 SAP

Spec 456065_R00 page 2/2

Printed Copies will not systematically be updated. Valid at the date of issue. Information is subject to change.