

Kremsmünster, 09.03.2020

VACUETTE® Urine Collection System – Product Data Specifications

Item N°¹	724310
Description¹	<i>Urine Beaker with Integrated Transfer Device, 100ml</i>
Materials¹	<i>Beaker PP Lid PE Needle Stainless steel Sleeve Synthetic rubber</i>
Additive^{1,2}	<i>none</i>
Legal Manufacturer²	Greiner Bio-One GmbH Bad Haller Str. 32 A-4550 Kremsmünster
Sterilisation¹	sterile, SAL 10 ⁻³
Directive²	IVDD 98/79/EC
Shelf Life²	60 months
Packaging¹	12 bags of 50 pcs. 600 pcs. In total

Sources:

1 AT_VAC_SAMP_724310_UrineBeaker_rev02

2 SAP

Spec 724310_R00

Printed Copies will not systematically be updated. Valid at the date of issue. Information is subject to change.

page 1/2

Kremsmünster, 09.03.2020

Label information¹

Urine Beaker with Integrated Transfer Device **greiner bio-one**

100 ml, sterile
single-packed, 200 pcs. in total

REF 724310	CODE 128 (ALL CHARS)
LOT A0000000	CODE 128 (ALL CHARS)
yyyy-mm-dd	CODE 128 (ALL CHARS)

Prod. Ord. 0000000 CODE 128 (ALL CHARS) (Prod. order n°.) Cart. 0000 CODE 128 (ALL CHARS) (Cart. n°.)

(SAL 10⁻³)
 0000

GS1 UDI

GS1 UDI Code (ALL CHARS)

(01)00000000000000(17)000000(10)A000000K

Greiner Bio-One GmbH
Bad Haller Str. 32
4550 Kremsmünster
Austria

Production location:
Greiner Bio-One Hungary Kft.
Fertősor 7
9200 Mosonmagyaróvár
Hungary

HIBC
DATA-MATRIX

Made in Hungary

Urine Beaker with Integrated Transfer Device

REF 724310	GS1 UDI DATA-MATRIX (01)00000000000000 (17)000000 (10)A000000K
LOT A0000000	
yyyy-mm-dd	

Production location:
Greiner Bio-One Hungary Kft., Fertősor 7, 9200 Mosonmagyaróvár, Hungary

(SAL 10⁻³)
 Made in Hungary

Greiner Bio-One GmbH, Bad Haller Str. 32,
A-4550 Kremsmünster, Austria

Sources:
1 AT_VAC_SAMP_724310_UrineBeaker_rev02
2 SAP

Spec 724310_R00
Printed Copies will not systematically be updated. Valid at the date of issue. Information is subject to change.