

Kremsmünster, 17.01.2020

## VACUETTE® Urine Collection System – Product Data Specifications

<b>Item N°1</b>	450251
<b>Characterisation<sup>1,2</sup></b>	Urine Transfer Device, short
<b>Description<sup>1</sup></b>	<i>Dimension total length 12.5cm</i>
<b>Materials<sup>1,3</sup></b>	<i>Holder Polypropylene Suction tube Polyethylene Needle Stainless steel</i>
<b>Additive<sup>1,2</sup></b>	<i>none</i>
<b>Manufacturer<sup>2</sup></b>	Greiner Bio-One GmbH Bad Haller Str. 32 A-4550 Kremsmünster
<b>Sterilisation<sup>1</sup></b>	none
<b>Directive<sup>2</sup></b>	IVDD 98/79/EC “Directive 98/79/EC of the European Parliament and council of 27th of October 1998 on in vitro diagnostic medical devices”. Classification: Other device (all devices except Annex II and self-testing devices).
<b>Storage<sup>1</sup></b>	none
<b>Shelf Life<sup>2</sup></b>	none
<b>Packaging<sup>1</sup></b>	12 bags of 50 pcs. 600 pcs. In total

Sources:

1 AT\_SUP\_SAMP\_450251\_450252\_Urine Transfer Device\_Carton Label\_rev02

2 SAP

3 A\_AN\_SOP\_040302\_223\_Rev01


Spec 450251\_R00 page 1/2


Printed Copies will not systematically be updated. Valid at the date of issue. Information is subject to change.

Kremsmünster, 17.01.2020

### Label information<sup>1</sup>

### Carton – Examples




**Urine Transfer Device Short**  greiner bio-one

12 bags of 50 pcs., 600 pcs. in total


REF **450251** CODE 128 (ALL CHARS)

LOT **A0000000** CODE 128 (ALL CHARS)

Prod. Ord. CODE 128 (ALL CHARS) (Prod. order n°) Cart. CODE 128 (ALL CHARS) (Cart. n°) 00000

CE 

GS1 UDI  
GS1 UDI Code (ALL CHARS)  
(01)0000000000000(10)A000000K

 Greiner Bio-One GmbH  
Bad Haller Str. 32  
4550 Kremsmünster  
Austria

Production location:  
Greiner Bio-One Hungary Kft.  
Fertőszőlő 7  
6200 Mosonmagyaróvár  
Hungary

DATA-MATRIX  
HIBC  
Made in Hungary

Sources:

1 AT\_SUP\_SAMP\_450251\_450252\_Urine Transfer Device\_Carton Label\_rev02

2 SAP

3 A\_AN\_SOP\_040302\_223\_Rev01

Spec 450251\_R00 page 2/2

Printed Copies will not systematically be updated. Valid at the date of issue. Information is subject to change.