

Kremsmünster, 17.01.2020

VACUETTE® Urine Collection System – Product Data Specifications

Item N°1	450252
Characterisation^{1,2}	Urine Transfer Device, long
Description¹	<i>Dimension total length 20cm</i>
Materials^{1,3}	<i>Holder Polypropylene Suction tube Polyethylene Needle Stainless steel</i>
Additive^{1,2}	<i>none</i>
Manufacturer²	Greiner Bio-One GmbH Bad Haller Str. 32 A-4550 Kremsmünster
Sterilisation¹	none
Directive²	IVDD 98/79/EC “Directive 98/79/EC of the European Parliament and council of 27th of October 1998 on in vitro diagnostic medical devices”. Classification: Other device (all devices except Annex II and self-testing devices).
Storage¹	none
Shelf Life²	none
Packaging¹	12 bags of 50 pcs. 600 pcs. In total

Sources:

1 AT_SUP_SAMP_450251_450252_Urine Transfer Device_Carton Label_rev02

2 SAP

3 A_AN_SOP_040302_223_Rev01





Spec 450252_R00 page 1/2

Printed Copies will not systematically be updated. Valid at the date of issue. Information is subject to change.

Kremsmünster, 17.01.2020

Label information¹

Carton – Examples

Urine Transfer Device 	
Short greiner bio-one	
12 bags of 50 pcs., 600 pcs. in total	
REF 450251	CODE 128 (ALL CHARS)
LOT A0000000	CODE 128 (ALL CHARS)
Prod. Ord. <small>CODE 128 (ALL CHARS) (Prod. order n°)</small>	Cart. <small>CODE 128 (ALL CHARS) (Cart. n°)</small> 00000
 	
GS1 UDI	
GS1 UDI Code (ALL CHARS)	
(01)0000000000000(10)A000000K	
 Greiner Bio-One GmbH Bad Haller Str. 32 4550 Kremsmünster Austria	Production location: Greiner Bio-One Hungary Kft. Fertősziget 7 6200 Mosonmagyaróvár Hungary
HIBC DATA-MATRIX	
Made in Hungary	

Sources:

1 AT_SUP_SAMP_450251_450252_Urine Transfer Device_Carton Label_rev02

2 SAP

3 A_AN_SOP_040302_223_Rev01

Spec 450252_R00 page 2/2

Printed Copies will not systematically be updated. Valid at the date of issue. Information is subject to change.