



## VACUETTE® Blood Culture Holder + Luer Adapter



- ↪ Single-packaged and sterile
- ↪ Can be used with various size blood culture bottles and also with blood collection tubes
- ↪ Manufactured with a normed male Luer adapter
- ↪ No separate insert needed
- ↪ Latex-free, non-DEHP product
- ↪ Convenient pouch packaging, keeping waste to a minimum

# VACUETTE® Blood Culture Holder + Luer Adapter

Few collection errors are as costly to the hospital or the laboratory as blood cultures that are falsified during specimen-collection practices.

False positive results can be caused, for example by, contamination from the patient's skin, the equipment used to take the sample and transfer it to the culture bottle, the hands of the person taking the blood sample or the general environment.

By choosing a single packaged, sterile product, this risk is reduced to a minimum. With this new sterile, single-packaged **VACUETTE®** Blood Culture Holder + Luer Adapter we are helping medical staff to optimize blood culture collection and further reduce the risk of false positive results, thereby helping health systems around the world to save costs, and allowing patients to receive the best possible treatment and be spared unnecessary treatments and delays.

## Intended Use:

The **VACUETTE®** Blood Culture Holder + Luer Adapter is used for a connection with a female Luer system in order to collect venous blood in blood culture bottles and / or evacuated blood collection tubes. This device is to be used by properly trained healthcare professionals only in accordance with these instructions.



## Beneficial features of the product

- Single-packaged and sterile
- Can be used with various size blood culture bottles (wide neck or narrow neck, suitable i.e. for BacT/Alert® and BD BACTEC™ bottles) as well as with blood collection tubes
- Manufactured with a normed male Luer adapter (applying ISO 594), thereby fitting with any other normed female Luer connector
- No need for a separate insert
- Latex-free, non-DEHP product
- Convenient pouch packaging, keeping waste to a minimum

## Instructions for Use

For the complete handling recommendations please refer to the Instructions for Use at [www.gbo.com](http://www.gbo.com).



**1** Connect the device to a female Luer connector.



**2** Hold the blood culture holder firmly and press it down on the bottle.



**3** After obtaining the correct amount of blood, continue (if required) with blood collection into **VACUETTE®** Blood Collection Tubes.



**4** Dispose of in an approved sharps disposal container.

Item N°	Description	Packaging	
		Inner	Outer
450197	<b>VACUETTE®</b> Blood Culture Holder + Luer Adapter 20G, single-packed, sterile, latex-free	24 pcs.	240 pcs.