

Kremsmünster, 09.03.2020

VACUETTE® Blood Collection System – Product Data Specifications

| | | |
|---------------------------------------|---|----------------------------------|
| Item N°¹ | 450070 | |
| Characterisation^{1,2} | Luer Adapter | |
| Description¹ | <i>Gauge</i> | <i>20G</i> |
| Materials^{1,3} | <i>Container, Hub</i> | <i>PP</i> |
| | <i>Needle</i> | <i>Stainless Steel, Silicone</i> |
| | <i>Sleeve</i> | <i>Synthetic Rubber</i> |
| Additive^{1,2} | <i>none</i> | |
| Legal Manufacturer² | Greiner Bio-One GmbH Bad Haller Str. 32 A-4550 Kremsmünster | |
| Sterilisation¹ | Sterile SAL 10 ⁻⁶ | |
| Directive² | Classification according to MDD 93/42/EEC | |
| Storage¹ | 4-36°C | |
| Shelf Life² | 60 months | |
| Packaging¹ | 20 boxes of 100 pcs. → 2000 pcs. 3 cartons of 2000 pcs. → 6000 pcs. in total | |

Sources:

1 AT_SUP_SAMP_450070_Lueradapter_Carton_label_rev03, AT_SUP_SAMP_450070_LuerAdapter_InnerBox_Label_rev00

2 SAP

3 STED-60579-03-01 Rev.00

Spec 450070_R00

Printed Copies will not systematically be updated. Valid at the date of issue. Information is subject to change.

page 1/2

Kremsmünster, 09.03.2020

Label information¹

Carton Labels:

VACUETTE® greiner bio-one

Luer Adapter

20G, sterile, not made with natural rubber latex
20 boxes of 100 pcs., 2000 pcs. in total

REF **450070** CODE 128 (ALL CHARS)

LOT **A0000000** CODE 128 (ALL CHARS)

yyyy-mm-dd CODE 128 (ALL CHARS)

CE 0123 STERILE EO

Production location:
Hiro Medical Industries Ltd., Takabayashi Plant, 2-19-64, Matsubara, Takabayashi-eki, Gunma, 374-8510 Japan
Made in Japan, distributed by Greiner Bio-One GmbH, Austria

4°C 40°F 36°C 97°F

GS1 UDI

GS1 UDI Code (ALL CHARS)

(01)00000000000000(17)000000(10)A000000K

Greiner Bio-One GmbH
Bad Haller Str. 32
4550 Kremsmünster
Austria

HIBC DATA-MATRIX

CAUTION: U. S. Law restricts this device to sale by or on the order of a physician.
Medical device for venous blood collection and in vitro diagnostic use.

VACUETTE® greiner bio-one

Luer Adapter

20G, sterile, not made with natural rubber latex
3 cartons of 2000 pcs., 6000 pcs. in total

REF **450070** CODE 128 (ALL CHARS)

LOT **A0000000** CODE 128 (ALL CHARS)

yyyy-mm-dd CODE 128 (ALL CHARS)

CE 0123 STERILE EO

Production location:
Hiro Medical Industries Ltd., Takabayashi Plant, 2-19-64, Matsubara, Takabayashi-eki, Gunma, 374-8510 Japan
Made in Japan, distributed by Greiner Bio-One GmbH, Austria

4°C 40°F 36°C 97°F

GS1 UDI

GS1 UDI Code (ALL CHARS)

(01)00000000000000(17)000000(10)A000000K

Greiner Bio-One GmbH
Bad Haller Str. 32
4550 Kremsmünster
Austria

HIBC DATA-MATRIX

CAUTION: U. S. Law restricts this device to sale by or on the order of a physician.
Medical device for venous blood collection and in vitro diagnostic use.

Inner Box Label:

REF **450070**

LOT **A0000000** GS1 UDI

yyyy-mm-dd

HIBC DATA-MATRIX

(01)00000000000000(17)000000(10)A000000K

REF **450070**

LOT **A0000000**

yyyy-mm-dd

GS1 UDI

GS1 UDI Code (ALL CHARS)

(01)00000000000000(17)000000(10)A000000K

HIBC DATA-MATRIX

Sources:

1 AT_SUP_SAMP_450070_Lueradapter_Carton_label_rev03, AT_SUP_SAMP_450070_LuerAdapter_InnerBox_Label_rev00

2 SAP

3 STED-60579-03-01 Rev.00

Spec 450070_R00

Printed Copies will not systematically be updated. Valid at the date of issue. Information is subject to change.