




TECHNICAL DATA SHEET

product name	color	case quantity	case weight(kg)	case volume(m ³)
JS-TQ003	sky blue	1,000	7.50	0.021



SPECIFIC USE - DESCRIPTION

Product intended for the retention of blood flow by oppression of a body limb. To perform blood sample extractions, differentiate a vein, etc.

Manufactured with TPE, thermoplastic elastomer.

Product free of Latex, non-cytotoxic or irritants



Ningbo Jenius Polymer Sci-Tech Co.,Ltd.



Manufactured with TPE, thermoplastic elastomer.

Product free of Latex, non-cytotoxic or irritant.

Length: 450 ± 5 mm

Width: 25 ± 1 mm

Thickness: 0.635 ± 0.02 mm


Tension strength ≥ 3 MPa

Elongation $> 300\%$



TECHNICAL DATA SHEET

product name	color	case quantity	case weight(kg)	case volume(m ³)
JS-TQ002	sky blue	1,000	6.30	0.018



SPECIFIC USE - DESCRIPTION

Product intended for the retention of blood flow by oppression of a body limb. To perform blood sample extractions, differentiate a vein, etc.

Manufactured with TPE, thermoplastic elastomer.

Product free of Latex, non-cytotoxic or irritants

Manufactured with TPE, thermoplastic elastomer.

Product free of Latex, non-cytotoxic or irritant.


Length: 450 ± 5 mm

Width: 25 ± 1 mm

Thickness: 0.635 ± 0.02 mm

Tension strength ≥ 3 MPa

Elongation $> 300\%$

 Marketing Department	 Technical, Development and Quality Director	 Sales & Marketing Director
---	---	--