



# IMPROVE YOUR LAB WORK FLOW EFFICIENCY

## **VACUETTE®** UNICAP Decapper

The VACUETTE® UNICAP Decapper quickly and safely decaps VACUETTE® blood collection tubes. Tubes just need to be inserted into a "belt" with adapters of the same or different composition. No time-consuming sorting of tubes before decapping is required, since the UNICAP Decapper can open different types of VACUETTE® tubes in a single operation.

The VACUETTE® UNICAP Decapping Device can decap up to 40 VACUETTE® tubes of the dimensions 13 x 75, 13 x 100 and 16 x 100 mm in one run. The belt can be individually composed of 75 and 100 mm tube length adapter elements.

- / Quick and safe opening of VACUETTE® Tubes
- / Directly after inserting the tubes into an adapter chain
- / No presorting of the tubes necessary
- / Robust and reliable
- / Easy to use

### VACUETTE® UNICAP DECAPPER ACCESSORIES.

Greiner Bio-One offers two different belts, each consisting of 40 individual tube adapters:

One for tube length 75 mm and one for 100 mm.



## PRODUCT OVERVIEW

FOR FURTHER INFORMATION AND/OR  
A DEMO PLEASE CONTACT US.

### VACUETTE® UNICAP Decapper

*with two belt adapters (for 75 mm and 100 mm tubes)*

Item no.	Qty. inner / outer
855904	- /1

### VACUETTE® UNICAP Decapper Accessories

Item no.	Description	Qty. inner / outer
855907	Belt Adapter for VACUETTE® UNICAP Decapper for 100 mm tubes	- /1
855908	Belt Adapter for VACUETTE® UNICAP Decapper for 75 mm tubes	- /1