



www.Hettich.dk

DATA SHEET

Basic data

Product:	Sharps Container 2 Lt P-SERIES
Item code:	DP.10.124

Characteristics data

Nominal volume:	2 Lt
Usefull volume:	1,6 Lt
Body color:	Yellow
Lid color:	Yellow
Print color:	Black
Dimensions: (assembled)	143x143x174,5 mm
Top width:	143 mm
Bot diameter:	114 mm
Height:	174,5 mm
Aperture diameter:	72,6 mm
Weight:	144 g
Material:	Polypropylene
Included PCR:	Minimum 40%
Diposal:	To be incinerated
Temporary lock:	Yes
Permanent lock:	Yes
Handle:	Yes
Personalized label option:	Yes



Logistic data

Carton quantity:	72 pcs
Carton weight:	11,32 kg
Carton dimensions:	600x400x430 mm
Carton volume:	0,1032 m ³
Pallette quantity:	1 440 pcs
Pallette weight:	244 kg (including pallette)
Pallette dimensions:	120x80x235 cm

Certifications data

UN 3291:	Yes
ISO 23907-1:	DQS GC Mark ISO 23907-1:2019
Quality management system:	ISO 9001 ISO 14001 ISO 45001 ISO 13485
CE Mark:	Marking not required as product is not a medical device (EU directive 93/42/EEC)

