



www.Hettich.dk

DATA SHEET

Basic data

Product:	Sharps Container 5 lt R-SERIES
Item code:	DP.10.50
Manufacturer:	Dunavplast Korp d.o.o.

Characteristics data

Nominal volume:	5 lt
Usefull volume:	4,25 lt
Body color:	Yellow
Lid color:	Yellow
Print color:	Black
Dimensions: (assembled)	190x190x241 mm
Top width:	190 mm
Bot diameter:	162 mm
Height:	241 mm
Weight:	342 g
Material:	Polypropylene
Included PCR:	Minimum 40%
Diposal:	To be incinerated
Temporary lock:	Yes
Permanent lock:	Yes
Handle:	Yes
Personalized label option:	Yes



Logistic data

Carton quantity:	30 pcs
Carton weight:	11,40 kg
Carton dimensions:	600x400x430 mm
Carton volume:	0,1032 m ³
Pallette quantity:	600 pcs
Pallette weight:	243 kg
Pallette dimensions:	120x80x235 cm

Certifications data

UN 3291:	Yes
ISO 23907-1:	DQS GC Mark ISO 23907-1:2019
	ISO 9001
Quality management system:	ISO 14001
	ISO 45001
	ISO 13485
CE Mark:	Marking not required as product is not a medical device (EU directive 93/42/EEC)