



www.Hettich.dk

DATA SHEET

Basic data

Product:	Sharps Container 6 lt XL P-SERIES
Item code:	DP.10.173
Manufacturer:	Dunavplast Korp d.o.o.

Characteristics data

Nominal volume:	6 lt
Usefull volume:	5,1 lt
Usefull height:	up to 40 cm
Body color:	Yellow
Lid color:	Yellow
Print color:	Black
Dimensions: (assembled)	143x143x471,6 mm
Top width:	143 mm
Bot diameter:	114 mm
Height:	471,6 mm
Aperture diameter:	72,6 mm
Weight:	435 g
Material:	Polypropylene
Included PCR:	Minimum 40%
Diposal:	To be incinerated
Temporary lock:	Yes
Permanent lock:	Yes
Handle:	Yes
Personalized label option:	Yes



Logistic data

Carton quantity:	16 pcs
Carton weight:	7,65 kg
Carton dimensions:	400x400x550 mm
Carton volume:	0,0088 m ³
Pallete quantity:	288 pcs
Pallete weight:	153 kg
Pallete dimensions:	120x80x180 cm

Certifications data

UN 3291:	Yes
ISO 23907-1:	DQS GC Mark ISO 23907-1:2019
Quality management system:	ISO 9001 ISO 14001 ISO 45001
CE Mark:	Marking not required as product is not a medical device (EU directive 93/42/EEC)