



## DOSESAVER™ LD

### 1 ML SYRINGE WITHOUT NEEDLE

LOW DEAD VOLUME OF ABOUT 0,062 ML



# Dosesaver™ LD

## 1 ml syringe without needle

LOW DEAD VOLUME OF ABOUT 0,062 ML

- Available in 1 ml volume
- Grey gasket instead of the regular black gasket for better visibility of the scale
- Reduces risk of over-sampling
- Reduces drug waste

CE<sub>0123</sub>  
Class I  
Medical Devices  
Rule 1  
MDD 93/42/EEC  
UMDNS 13929

### Product Compliance

- Complies with the following norms, directives, and regulations:
  - ISO 7864:2016
  - EN ISO 80369-7:2021
  - EN ISO 11607-1:2019
  - EN ISO 11607-2:2019
  - ISO 80369-7:2021
  - ISO 7886-1: 2017
  - ISO 80369-7: 2021
  - ISO 11607-1: 2019
  - ISO 11607-2: 2019
  - ISO 10993-7: 2008/Amd 1: 2019

### Manufacturing Details

**Legal manufacturer:** Nipro Corporation

**Country of origin:** Thailand

### Sterilization and Shelf Life

#### Sterilization

EtO (Ethylene oxide)



#### Single-use only



#### Shelf life

5 years



### Materials Used



# PRODUCT RANGE OVERVIEW

Reference code	Description
SY3-LD-1SC-ET	LD 3 part syr 1 ml W/O needle slip

**Not made with latex, DEHP, PVC or BPA**

## Packaging Details

Individual unit			Inner box				Outer box				
Length (mm)	Width (mm)	Weight (g)	Length (mm)	Width (mm)	Height (mm)	Quantity	Length (mm)	Width (mm)	Height (mm)	Weight (kg)	Quantity
140	96	4,62	202	135	91	100	229	482	337	4,62	1000

## Languages

27 languages on all levels of packaging:

Arabic (AR), Bulgarian (BG), Czech (CS), Danish (DA), English (EN), German (DE), Greek (EL), Spanish (ES), Estonian (ET), Finnish (FI), French (FR), Croatian (HR), Hungarian (HU), Italian (IT), Lithuanian (LT), Latvian (LV), Dutch (NL), Norwegian (NO), Polish (PL), Portuguese (PT), Romanian (RO), Russian (RU), Slovakian (SK), Slovenian (SL), Serbian (SR), Swedish (SV), Turkish (TR)



## Storage conditions






Do not store at extreme temperature and humidity.

## Transport conditions

Closed and dry.

## Label details

	Individual unit	Inner box	Outer box
Product name	✓	✓	✓
Product size	✓	✓	✓
Quantity (pieces)	—	✓	✓
REF	—	✓	✓
LOT	✓	✓	✓
	✓	✓	✓
STERILE EO	✓	✓	✓
	✓	✓	✓

	Individual unit	Inner box	Outer box
CE	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
EC REP	✓	✓	✓
	—	✓	✓
Made in Thailand	✓	✓	✓

